

Technical Data

Part No.: 501592V07001000
Product: Turbomolecular Pump System TURBOLAB 450
consisting of:
Turbomolecular Pump TURBOVAC 450 i
Diaphragm Vacuum Pump DIVAC 3.0 T

Inlet connection:	DN 160 CF
Outlet connection:	Silencer, G 1/8" Inside thread
Pumping speed	
N ₂ - Nitrogen:	430 l/s
Ultimate pressure:	1.0 x 10 ⁻¹⁰ mbar
Mains connection:	1-ph, 110 V ±10%, 50/60 Hz 1-ph, 230 V ±10%, 50/60 Hz
Max. Power consumption:	240 W
Dimensions:	see dimension sheet
Weight:	≈ 35.0 kg

Technical data are subject to change